



CORONAVIRUS PREVENTION



**SANITISING AND
SAFETY PRODUCTS
FOR COVID-19**

PPE KIT



- Fabric used for Coverall is SITRA Certified.
- Protective coverall is of 85GSM Non woven laminated
- Includes N95 Mask, Hood cover, Goggles, Gloves, Shoe cover, safety disposal bag.
- Sterilization process is conducted appropriately
- Daily capacity of 3000 PPE kit and can be increased in case of bulk quantity order
- ISO , CE Certificate and Drug license available.



Vivo V15
AI Triple Camera

PPE KIT

Clinisafe

Personal Protection kit.

ISO 16603:2004

Coverall Resistance to Penetration of Blood and body fluids using synthetic blood.

Mask 3ply



Goggles

Wear on optical



Coverall Gown

Hood cover

Leggings

**DRDO OFB
Approved coverall**

**White / Blue colour
70+ Gsm lamination gown
With seam sealing tapes**

Nitrile Gloves



**Bio Hazard Bag for
Waste collection**



Product By:

- DRDO and ISO 1660:2004 Certified PPE Kit.
- Protective coverall consist of 5 layer , consist of 2 layer of Spunbond polypropylene, 2 layer of melt down polypropylene, 1 layer of spun bond.
- Includes 3ply Mask, Hood cover, Goggles, Gloves, Shoe cover, Safety disposal bag
- Fabric used is 70gsm SSMMS which is heavy quality then Non-woven 80GSM.
- Sterilization process is conducted appropriately.
- Daily capacity of 1000 PPE kit and can be increased in case of bulk quantity order.

N95 Mask



Total Five Layer:

- Outer Most – 70 GSM Non Woven
- One 25 GSM Non Woven
- And Three 20 GSM 98.9% BFE

Compliant melt blown inside the mask.

Certificate

Application for
Personal Protective Equipment (PPE) Directive 89/686/EEC
This is to certify that the Quality System for the Product Group:

Manufacturing and Supply of N95, KN95 and 3PLY Mask with Meltblown and Spunbond Material, Nonwoven Full Body Cover and Shoes Cover for PPE Kit

Defined by the manufacturer as per the PPE
A Product of

Complies with the requirement of the
"PPE"
The quality system for these products has been assessed according
to the procedure of conformity assessment Prescribed.
This certificate is issued in pursuant to the assessment done,
under the provision

Personal Protective Equipment Directive 89/686/EEC
QCC has liabilities limited to the audit for certification of the products and in case of
alteration of the product, not agreed upon by us, this declaration will lose its validity.

Certificate No. : CE/016014/0320
Original Certificate Date : 23-March-2020
Issue Date : 23-March-2020
Expiry Date : 22-March-2023
To check this certificate status visit:
"http://uasi.uk.com/certifiedorganization.html"

Authorized Signature
For Quality Control Certification
UK Office: 1929, Chynoweth House,
Trevelyan Park, Truro (TR40LN, Cornwall, UK
India Office: 2nd Floor, Anand Market,
Narada Mandi, Delhi - 110 040, India

Quality Control Certification awarded by UASI, UK
This certificate doesn't provide the certified organization with immunity from its legal obligations.
The certificate is issued on the basis of the information provided by the manufacturer.

LILABA ANALYTICAL LABORATORIES
COMPREHENSIVE ANALYTICAL SERVICES

Report No : LAL-240420-1

Issued To : [Redacted] Report Date : 25/04/2020
Lab ID No. : 240420-1
Your Ref. : [Redacted]

PART A : PARTICULARS OF SAMPLES

1. Date Of Receipt : 24/04/2020	12. Declared Value if any : N/A
2. Analysis Starting Date : 24/04/2020	13. Mode of Packing : Individual
3. Analysis Completion date : 25/04/2020	14. Grade/Variety/Size : [Redacted]
4. Batch No. : [Redacted]	15. Sample Drawn By : Customer
5. Date of Mfg. : [Redacted]	16. Brand Name, if any : N/A
6. Date Of Expiry : [Redacted]	17. BSE Seal : [Redacted]
7. Test Report As Per : EN 14683	18. IO's Signature : N/A
8. Description of the Sample : [Redacted]	
9. Sample Name : N95 Face Mask	
10. Sample Qty. Received : 10	
11. Any other Information : [Redacted]	

PART B : TEST RESULT

Sr. No.	TEST NAME	SPECIFIED REQUIREMENT	RESULTS
1	Type	5 Layer	5 Layer
2	Color	White	White
3	Check for Integrity	Proper Ultrasonically Welded Layer	Complies
4	Absorbance Test	NMT 0.3 (000)	Complies
5	Residue on Evaporation	NMT 5 mg	Complies
6	Acidity or Alkalinity	NMT 1 ml of 0.001 M HCl/NaOH	Complies
7	Heavy metals ions Test	NMT 1 ppm	Complies
8	Reducing (Oxidizable) Test	NMT 2 ml of 0.002 M KMnO4	Complies
9	Breathing Resistance	NMT 40.0 l/min/cm2	Complies
10	Bacterial Filter Efficiency (BFE)	95.2%	98.2%
11	Particle Filter Efficiency (PFE)	NLT 95.2%	98.2%
12	Splash Resistance (mm Hg)	NLT 12.0	Complies
13	Toxicity	Non Toxic	Complies

Checked By : [Signature] Authorized Sign : [Signature]

Address : 2nd Floor, Galaxy Point, Above Hotel Anand, Opp. Anand Zoo, Sarthana Jakaraka, Varanasi Road, Surat-395006. (Dist.) India. Admin Contact : +91 78785 22920 / +91 97377 04415 Sales Inquiry : +91 96243 22956 / +91 97375 44666 Email: info@lilaba@gmail.com

LILABA ANALYTICAL LABORATORIES
COMPREHENSIVE ANALYTICAL SERVICES

Report No : LAL-240420-1

Remarks :
Sample confirm as per EN 14683 with respect to above tested Parameters only.

NOTE :
Please Refer Overleaf

Face Shield



← Anti-stress foam

350 micron
clear film →



American National Standards Institute
ANSI STANDARD

QUALITY CERTIFICATION



Certificate of Registration

This is to certify that The Quality Management System of Medical Devices



has been assessed and found to be in compliance with the requirements of the standard

ISO 13485:2016

for the following scope :

Manufacturing and Supply of Face Shield, Face Protector, PP Box, Pet Boxes, Safety Goggles, Disposable Face Shield for P.P.E. Kit

CERTIFICATE No. : 20ZGA04220M

ISSUED DATE : 25/03/2020
EXPIRY DATE : 24/03/2023

1ST SURVEILLANCE : 24/03/2021
2ND SURVEILLANCE : 24/03/2022



Authorised Signatory

INTERNATIONAL QUALITY CERTIFICATION SERVICES UK LTD
272, Bath Street, Glasgow, G2 4JR, U.K.

This Certificate is Intellectual Property of IQCS and can be maintained through surveillance and renewal audits.
Certificate should be returned to IQCS in case of non compliance of certification procedures.
Authority of this certificate can be verified at www.iqcert.co.uk / www.iqcert.co.uk
The Registration is not a Product Quality Certificate.



Nitrile Gloves



- a) Daily One time use Nitrile Gloves.
- b) ISO Certified , Latex Free.

Disinfectant Sanitizer Spray Machine.



Sodium Hypochlorite

ABOUT THE PRODUCT

SODIUM HYPOCHLORITE, COMMONLY KNOWN AS BLEACH, IS MOST FREQUENTLY USED AS A DISINFECTING AGENT. IT IS A BROAD-SPECTRUM DISINFECTANT THAT IS EFFECTIVE FOR THE DISINFECTION OF VIRUSES, BACTERIA, FUNGI, AND MYCOBACTERIUM.

SPECIFICATION

APPEARANCE : COLOURLESS TO PALE YELLOW LIQUID
ACTIVE CHLORINE CONTENT : 5-10%
FREE ALKALI (AS NaOH) (G/L) : MIN 5 (FOR 5% SOLUTION)
MIN 10 FOR (10% SOLUTION)

APPLICATION

IN VIEW OF THE ONGOING COVID-19 VIRUS, GOI HAS RECOMMENDED THE USAGE OF DILUTE SODIUM HYPOCHLORITE FOR DISINFECTING PUBLIC PLACES LIKE BUSES, HOSPITALS, RESTAURANTS, LOCAL STREETS, OFFICES, CONSTRUCTION SITES ETC.

IT IS ALSO USED AS A POWERFUL BLEACHING AGENT IN THE PULP & PAPER AND TEXTILE FIBRES.

ANOTHER IMPORTANT PROPERTY OF SODIUM HYPOCHLORITE IS THAT IT IS USED AS AN OXIDIZING AGENT, A PROPERTY MUCH EXPLOITED IN VARIOUS CHEMICAL REACTIONS.



50 KG PACKING OF SODIUM HYPOCHLORITE

DOSE (recommended)

FOR OPEN SPACES, DILUTE :-

-10 LITRE 5% SODIUM HYPOCHLORITE SOLUTION IN 50 LITRE WATER.

-10 LITRE 10% SODIUM HYPOCHLORITE SOLUTION IN 100 LITRE WATER.

NOTE:-

AS PER THE GUIDLINES OF THE GOI, AJANTA GROUP DOESNT RECOMMEND USING SODIUM HYPOCHLORITE ON "HUMANS".

PACKAGING

- M.S.R.L ROAD TANKERS

SPECIAL SMALLER PACKAGING CAN ALSO BE CONSIDERED ON REQUEST

Sanitizer E-Dispenser

E-SENSE

Auto Touchless Soap Dispenser



Place your hand
under the sensor



Advantages

- Prevents cross infection.
- Built-in infrared smart sensor.
- Fully automatic and touch-free operation.
- Eliminates countertop mess and wastage.
- Convenient 5 adjustable levels
- Ideal for lotions, liquid soap or sanitizer etc.

Pedal Sanitizer Stand

